

RAPID | SIMPLE | SELF-CONTAINED | LOW-COST | RELIABLE

A Powerful MDx Platform for Nucleic Acid and Protein Molecular Testing

OVERVIEW

For over 5 years, Akonni Biosystems has been a pioneer in the MDx industry, helping evolve the field of low-cost microarray technology for applications in clinical and life science research markets. Providing systems for technicians and researchers in academia, government, and industry labs and addressing a wide range of applications, Akonni Biosystems is proud to provide quality instrumentation, software solutions, and responsive technical support.

INSTRUMENTS AND SOFTWARE AVAILABLE

TruDiagnosis® System: Akonni TruDiagnosis Systems provide a workflow-integrated platform for multiplexed MDx and immunoassays. The system is modular in design and has a variety of configurations to meet the broad range of testing demands in most laboratory environments.

100-10000 Akonni TruDiagnosis System for multiplexed molecular testing. Includes Akonni TruDx™ 1000 Reader, Akonni TruCycler™ Satellite Thermal Cycler, and laptop computer with data collection and analysis software. Systems include one year warranty against defective parts. Accepts international AC; 90-264V, 50-60 Hz AC.

TruDx 1000 Reader: Akonni TruDx 1000 Readers provide a user-friendly image capture and analysis system optimized to read TruArray® tests. The reader provides exceptional value in a truly high performance yet portable unit.

100-10010 TruDx 1000 Reader for TruArray image capture and data analysis on a laptop computer. Readers include one year warranty against defective parts. Accepts international AC; 90-264V, 50-60 Hz AC.



TruCycler Satellite Thermal Cycler: Akonni TruCycler Satellite Thermal Cyclers are high performance thermal cyclers with an add-on slide block for performing TruArray hybridizations.

100-10020 TruCycler Satellite Thermal Cycler Instrument for TruArray test amplification and hybridization. Instruments include one year warranty against defective parts. Accepts international AC; 90-264V, 50-60 Hz AC.

To order call +1 301.698.0101

info@akonni.com | sales@akonni.com | support@akonni.com

For Research Use Only. Not for use in diagnostic procedures. No claim or representation is intended for the diagnosis, prevention, or treatment of a disease.

Protected by 1 or more U.S. Pat: 6,460,094; 6,442,734; 6,795,949; 7,096,435. Other patents issued and pending.

TruDiagnosis and TruArray are registered trademarks and the Akonni logo, TruDx, and TruCycler are trademarks of Akonni Biosystems, Inc.

© 2009 Akonni Biosystems. All rights reserved.