

FAST | EASY | AFFORDABLE | SAFE | RELIABLE

## Starter Kits

### OVERVIEW

Nearly every life science research and molecular diagnostic laboratory requires a rapid and robust means of extracting RNA and/or DNA. Akonni TruTip Kits deliver inhibitor-free, PCR-ready nucleic acid in as few as 4 minutes. TruTips can process a wide range of sample sources (blood, culture, sputum, urine, etc.), sample viscosities (thick or thin), and types (bacterial, viral or genomic).

With TruTip Starter Kits, we make it easy and affordable to get started. TruTip Starter Kits provide enough TruTips and buffers to perform 60 extractions for 2-ml kits, and 96 extractions for 1-ml kits. The 2-ml format accommodates larger samples, while the 1-ml format is compatible with multi-channel pipettes. For more viscous samples, select TruTip LPT tips; for less viscous samples and higher recovery, select TruTip SPT tips. As an added bonus, with a TruTip Starter Kit, you receive a RAININ® 100–1000µL (for 1-ml kits) or 200–2000µL (for 2-ml kits) EDP3-Plus™ Advanced Electronic Pipette with LTS® and a wall power supply.

### STARTER KITS AVAILABLE

**TruTip DNA Kits:** To extract DNA from bacterial or viral pathogens.

- 300-10606 TruTip DNA Starter Kit (2-ml Rainin LTS, 30 SPT/30 LPT)
- 300-10616 TruTip DNA Starter Kit (1-ml Rainin LTS, 48 SPT/48 LPT)

**TruTip RNA Kits:** To extract RNA from bacterial or viral pathogens.

- 300-10806 TruTip RNA Starter Kit (2-ml Rainin LTS, 30 SPT/30 LPT)

**TruTip XT DNA Kits with Beads:** To extract DNA from bacterial or viral pathogens that are tougher to lyse.

- 300-10306 TruTip XT DNA Starter Kit (2-ml Rainin LTS, 30 SPT/30 LPT)
- 300-10316 TruTip XT DNA Starter Kit (1-ml Rainin LTS, 48 SPT/48 LPT)

**TruTip DNA /RNA Blood Kits:** To extract nucleic acid from blood-borne bacterial or viral pathogens in samples up to 0.1 ml.

- 300-10404 TruTip DNA/RNA Blood Starter Kit (2-ml Rainin LTS, 60 SPT)



**TruTip DNA Kits for MTB:** To extract DNA from *Mycobacterium tuberculosis* culture or sputum of up to 1.0 ml or greater.

- 300-10105 TruTip Starter Kit for MTB (2-ml Rainin LTS, 60LPT)

**TruTip gDNA Saliva Kit:** To extract genomic nucleic acid from microliter quantities of liquid saliva, though readily adaptable to larger sample input volumes.

- 300-20106 TruTip gDNA Saliva Starter Kit (2-ml Rainin LTS, 30 SPT/30 LPT)

**TruTip gDNA Blood and Semen Kit:** To extract genomic nucleic acid from several hundred microliters of liquid blood or neat semen samples, though readily adaptable to smaller or larger sample input volumes.

- 300-20206 TruTip gDNA Blood/Semen Starter Kit (2-ml Rainin LTS, 30 SPT/30 LPT)

**Special Offer:** Get any Starter Kit starting at \$499\*, and receive a RAININ EDP3-Plus Advanced Electronic Pipette with LTS. Offer expires 6/30/2010.

To order call +1 301.698.0101

info@akonni.com | sales@akonni.com | support@akonni.com

\*Kits with 96 1-ml tips are \$699. For first time TruTip customers only, with a limit of one TruTip Starter Kit per account. Offer valid through June 30, 2010 and limited to US customers only.

For Research Use Only. Not for use in diagnostic procedures. No claim or representation is intended for the diagnosis, prevention, or treatment of a disease.

Patent pending: Akonni TruTip Extraction Tips.

The Akonni logo and TruTip are trademarks of Akonni Biosystems, Inc. RAININ and LTS are registered trademarks, and EDP3-Plus is a trademark of Rainin Instrument, LLC.

© 2009 Akonni Biosystems. All rights reserved.

TTP-DS-001 1004.2